



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LIKAMED NOCO X 330 treatment chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Emergency power supply: via 4 new built-in batteries 6 V / 1.2 Ah,
- Frequency: 50/60 Hz,
- **Medically certified upholstery in blue,**
- **The LIKAMED NOVO X 330** treatment chair provides free adjustment with the possibility of positioning in a wide range,
- The chair has excellent mobility, making it suitable for free movement around medical facilities,
- **Technical specifications:**
 - **Adjustment via wired handheld remote control:**
 - Seat height adjustment range: 57 – 82 cm,
 - Possibility of free adjustment of the backrest tilt within the range of 70° (**automatic synchronization of the leg rest with the tilt of the chair**) ,
 - Free adjustment of the leg rest inclination in the range of: 28°,
 - Programmed button for lying position,
 - Programmed button for sitting position,
 - Double-sided armrests with the possibility of manual free rotation,
 - Manually adjustable patient step extension in the range of 30 cm,
 - Mobile chair with built-in double-sided brake,
- **Remote control with built-in function lock (can be unlocked using the dedicated magnet included in the set),**
- **Included in the set:**
 - Wired remote control,
 - Magnet to unlock the remote control functions,
- **Dimensions:**
 - Headrest: 52 x 25 cm,
 - Backrest: 60 x 62 cm,
 - Seat: 62 x 50 cm,
 - Leg rest: 60 x 65 cm,
 - Total: 90 x 180 x 158 cm,
- Maximum load capacity: 175 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

