



DRAGER VAPOR 2000 inhalation agent vaporizer (Isoflurane) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product**,
- German production,
- **Drager Vapor 2000** is a safe and precise vaporizer for inhalants,
- **Product Features:**
 - The volume of the vaporizer is 300 ml, so the entire contents of a standard 250 ml anesthetic bottle can be poured into the vaporizer, even if there is anesthetic left inside from the previous refill,
 - The vaporizer does not require calibration during use,
 - Unique transport function - the vaporizers can be carried in any position, as the transport lock prevents the anesthetic from spilling, so it can be used immediately after transport,
 - All filling systems are designed specifically for specific inhalation anaesthetics, eliminating the risk of accidentally using the wrong anaesthetic agent.
 - Mounting system: The design of the vaporizer is very flexible. In addition to the Dräger DW-2000 connection system, the S-2000 mounting is also available, compatible with the Selectatec system. Compatible with the ISO 23 mm conical connection and the Dräger AutoExclusion mounting system,
 - **In this evaporator the SELECTATEC S-2000 mounting is used,**
 - Infusion type: Drager Fill,
 - Possibility to adjust the isoflurane concentration in the range of: 0.2 - 6%,
- **Included in the set:**
 - English user manual (PDF),
- Dimensions: 11 x 25 x 20 cm,
- Weight: 8 kg,
- Has a valid Technical Passport issued,
- Guarantee:
 - 6 months for the domestic market - Poland,
 - 3 months for international market,
- Financing options (Poland only): Installments, Leasing, Loan,

