



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER VAPOR 2000 inhalant vaporizer (Isoflurane/Isoflurane) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product

German manufacture,

Drager Vapor 2000

Is a safe and precise vaporizer of inhalants,

Product Features:

The volume of the vaporizer is 300 ml, so the entire contents of a standard 250 ml anesthetic bottle can be poured into the vaporizer, even when there is anesthetic left inside from the previous filling,

The vaporizer does not require calibration during use,

Unique transport function - vaporizers can be carried in any position, as the transport safety device prevents spillage of the anesthetic, so it can be used immediately after transport,

All filling systems are specially designed for specific inhalation anesthetics, which means eliminating the risk of mistakenly using the wrong anesthetic,

Attachment system: the design of the vaporizer is highly flexible. In addition to the Dräger DW - 2000 connection system, the S - 2000 attachment, compatible with the Selectatec system, is also available. Compatible with the ISO 23 mm conical connection and the Dräger

AutoExclusion attachment system,

The SELECTATEC S-2000 mount is used in the particular vaporizer,

Fill type: Drager Fill,

Adjustable isoflurane concentration in the range: 0.2 - 6 % ,

Included:

User manual in English (PDF),

Dimensions: 11 x 25 x 20 cm,

Weight: 8 kg,

Has a current Technical Passport issued,

Warranty:

6 months for the domestic market - Poland,

3 months for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

