



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### Medical assistant METRO FLEXLINE CRASH TROLLEY (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- Made in the USA,
- Mobile medical assistant METRO FLEXLINE CRASH TROLLEY with pull-out drawers,
- Technical specifications:
  - 4 pull-out drawers with the possibility of locking them,
  - Extendable IV stand with a total length of 101 cm,
  - Side bottle mount with gradual adjustment,
  - Mobile assistant on 4 braking wheels,
- Dimensions:
  - Top plate: 66 x 40 cm,
  - 3 pull-out drawers: 51 x 43 x 15 cm,
  - 1 pull-out drawer: 51 x 43 x 7 cm,
  - Total: 75 x 50 x 150 cm,
- Weight: 45 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
  - 6 months for domestic market (Poland),
  - 3 months to international market,
  - Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
- Financing options (Poland only): Installments, Leasing, Loan,

