



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEEC W100 Laboratory Incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- UK production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 240 W,
- **The LEEC W100** incubator is a spacious laboratory incubator with a capacity of 100 l.,
- Temperature range: up to 60°C
- Made of powder coated stainless steel,
- The easy-to-use temperature controller allows for quick programming of the device,
- The too low and too high temperature alarm is triggered when the set temperature is exceeded by 2° C,
- An independent, fail-safe switch protects the device from overheating,
- Possibility of setting the limit temperature by the operator,
- Air circulation is supported by a fan, which ensures even temperature distribution throughout the chamber,
- Temperature recovery within 2 minutes minimises heat loss when the door is opened,
- External doors equipped with a magnetic seal,
- Additional internal glass doors for the upper and lower parts of the chamber allow you to observe the process without the risk of the chamber suddenly cooling down,
- **Included in the set:**
 - Power cord,
 - 2 pcs. – trays 45 x 51 cm,
- Dimensions: 68 x 62 x 58 cm,
- Weight: 54 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

