



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT BE 400 laboratory incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 800 W,
- **The MEMMERT BE 400** incubator is ideal for applications requiring precise temperature control in the range from +20 °C to +70 °C,
- **The MEMMERT BE 400** with a capacity of 53 liters is designed for laboratories and applications requiring precise and reliable thermal control. The incubator is designed to provide exceptional performance in a wide range of applications, including cell culture, microbiological testing and environmental research.
- Thanks to natural convection and intelligent PID controller, it ensures even heat distribution,
- Stainless steel construction ensures easy maintenance and durability,
- **Technical specifications:**
 - Internal capacity: 53 l.,
 - Temperature range: from 20 °C to +70 °C,
 - Overheating protection: Triple safety system (switches off the heating at a temperature 3 °C higher than the set temperature, digitally adjustable temperature controller, mechanical limiter at a temperature 10 °C higher than the maximum),
 - Sensors: Pt100 class A sensor in a 4-wire circuit, mutual monitoring,
 - Temperature accuracy: 0.1 °C (for set and actual value),
 - Timer: up to 999 hours,
 - Heating elements: four-sided heating with a large surface area,
- **Included in the set:**
 - Power cord,
 - 2 pcs. – inner tray dimensions: 28 x 40 cm,
- Dimensions:
 - External: 55 x 48 x 67 cm,
 - Internal: 40 x 32 x 40 cm,
- Weight: 35 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

