



## MEMMERT S10 Laboratory Incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 900 W,
- **The MEMMERT S10** incubator is ideal for applications requiring precise temperature control in the range from +30 °C to +220 °C,
- **The 30-liter MEMMERT S10** is designed for laboratories and applications requiring precise and reliable thermal control. The incubator is designed to provide exceptional performance in a wide range of applications, including cell culture, microbiological testing and environmental research.
- Thanks to natural convection and intelligent PID controller, it ensures even heat distribution,
- Stainless steel construction ensures easy maintenance and durability,
- **Technical specifications:**
  - Internal capacity: 30 l.,
  - Timer: 120 minutes,
  - Temperature range: from 30 °C to +220 °C,
  - Heating elements: four-sided heating with a large surface area,
- **Included in the set:**
  - Power cord,
  - 3 pcs. – inner tray dimensions: 39 x 23 cm,
- Dimensions:
  - External: 51.5 x 44.5 x 56.6 cm,
  - Internal: 39.5 x 24.5 x 28.5 cm,
- Weight: 50 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for domestic market (Poland),**
  - **3 months to international market,**
  - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

