



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE VIO 300D surgical diathermy (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- **Software version:** 2.7.0,
- **Erbe Vio 300 D** is a surgical diathermy used in gastroenterology, endoscopy and pulmonology,
- **Features and benefits of the system:**
 - Modular structure of the system,
 - Color display,
 - Voltage regulation for economical, repeatable cutting and coagulation,
 - Adjustment of the electric arc - sparks, for cutting and coagulation with high energy and cutting in liquids,
 - Adjusting the output power to ensure a constant level during coagulation and devitalization,
 - NESSY – neutral electrode safety system,
 - Plug & play system – automatic instrument recognition,
 - Re-Mode function - remote control using a handle or foot attachment directly from the sterile operating field,
 - PreView function – preview function of the approximate effect on the tissue using the current settings,
 - The FocusView function reduces the displayed information to the most important ones, showing only the operating parameters of the instrument currently in use,
- **Monopolar socket:**
 - **Coagulation mode:**
 - SOFT COAG,
 - SWIFT COAG,
 - FORCED COAG,
 - COAG SPRAY,
 - CLASSOC COAGO,
 - **Cutting modes:**
 - AUTO CUT,
 - HIGH CUT,
 - DRY CUT,
 - ENDO CUT,
 - ENDO CUTQ,
 - **Bipolar socket:**
 - **Coagulation mode:**
 - BIPOLAR SOFT,
 - BIPOLAR SOFT (AUTO STOP),
 - BIPOLAR FORCED,
 - **Cutting mode:**
 - BIPOLAR CUT,
- **Specification:**
 - Maximum cutting power: 300 W at 500 Ω,
 - Maximum coagulation power: 200 W at 500 Ω,
 - Maximum bipolar power 370 W at 500 Ω,
- **Possibility to select menu language:** ENG, GER, NOR, SVE, DUN, GREK, NED, FR, IT, ESP, PT, RUS, PL, CS, JP, CHIN, TUR,
- **Included in the set:**
 - Brand new JackKNIFE BOWA monopolar electrode cable - 2 buttons, 4 mm shaft, , 4.5 m cable), (no electrodes included, can be retrofitted for an additional fee),
 - Brand new BOWA 4 m bipolar electrode cable (can be retrofitted with bipolar forceps for an additional fee),
 - Brand new cable with BOWA clip for disposable neutral electrodes (electrodes not included, can be retrofitted for an additional fee),
 - ERBE wired foot manipulator,
 - ERBE wired foot switch,
 - Mobile diathermy trolley **ERBE VIO-CART** ,
 - English user manual (PDF),
 - Power cord,
- Overall dimensions: 60 x 64 x 120 cm,
- Weight: 50 KG,
- Cat 02,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

