



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT WB10 water bath (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 1200 W,
- **Technical specifications:**
 - The bath is made entirely of high-quality acid-resistant steel,
 - Usable capacity: 10 liters,
 - Temperature control in the range: from ambient temperature to 95°,
 - Timer adjustment range: from 1 minute to 99 hours,
 - Built-in digital display,
 - Simple and intuitive operation,
- Overall dimensions: 35 x 21 x 14 cm,
- Weight: 14 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

