



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Cryo-T Elephant Mini Cryotherapy Apparatus (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production, (2024),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 510 W,
CRYO-T Elephant Mini

is an apparatus for cryotherapy with nitrogen, made with heater technology. It is an apparatus in a housing with a mobile base adapted to accommodate a 35-liter tank of liquid nitrogen,

The Elephant mini apparatus achieves low temperatures in a completely different way than apparatuses previously on the market. The new technology in this regard is based on the use of the heat of evaporation of liquid nitrogen. Thanks to modern design solutions, much more economical operation has been achieved,

The apparatus is controlled by a large touch panel with an LCD screen on which all the parameters of the treatment are shown. The intensity of cooling is achieved by gradual adjustment of the blowing force, selecting the optimal treatment parameters for the size of the body area subjected to cryostimulation,

Features:

Quick start:

Quick start allows to proceed immediately after starting the device. In a short period of time, the device produces cold between -110°C - 160°C, necessary for effective treatment. The target temperature is reached 1 minute after pressing the START button,

Patented technology:

State-of-the-art cryogenic tube guarantees unparalleled comfort and total freedom of movement (protected by international patent). It does not stiffen during the procedure, making the operator's hands fatigue-free. The insulating material used in the cryogenic cable prevents it from freezing,

Savings of up to 40% :

The Cryo-T® Elephant is equipped with innovative heater technology, guaranteeing significant reductions of up to 40% in refrigerant consumption ,

Non-freezing hose;

The Cryo-T Elephant's patented technology features a non-freezing, patented cryogenic hose for an unparalleled level of convenience in use and total freedom of movement,

Fast and safe treatment:

Treatment time in one area of the body ranges from 30 seconds to 5 minutes, while in the facial area it takes 2 to 3 minutes. If the patient has an increased amount of fat or muscle, the treatment duration can be extended to 10 minutes. Correct application of cold depends on the experience and knowledge of the doctor or therapist, as well as precise and continuous monitoring of the tissues' response to the low temperature,

Application:

Physiotherapy and rehabilitation,

Sports medicine,

Aesthetic medicine,

Wellness,

Rheumatology,

SPA,

Technical specifications:

Refrigerant: LN2,

Cold jet temperature: up to -160°C,

5-level jet intensity control,

Built-in coolant level indicator,

Included:

Power cord,

Protective gloves,

Safety glasses,

Overall dimensions: 55 x 65 x 150 cm,

Weight: 70 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

