



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEMMERT BE 500 Laboratory Warming Machine (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 900 W,  
Incubator  
MEMMERT BE 500  
is ideal for applications requiring precise temperature control from +30 °C to +60 °C,  
MEMMERT BE 500  
with a capacity of 108 liters, is designed for laboratories and applications requiring precise and reliable thermal control. The incubator is designed to provide exceptional performance in a wide range of applications, including cell culture, microbiological testing and environmental research.  
With natural convection and an intelligent PID controller, it provides even heat distribution,  
Stainless steel construction ensures easy maintenance and durability,  
Technical specifications:  
Internal capacity: 108 l,  
Temperature range: from 30 °C (min. 5 °C above ambient temperature) to +70 °C,  
Protection against overheating: Triple protection system (shuts off heating at 3 °C above the set temperature, digitally controlled temperature controller, mechanical limiter at 10 °C above the maximum temperature),  
Sensors: Pt100 class A sensor in 4-wire circuit, mutual monitoring,  
Temperature accuracy: 0.1 °C (for setpoint and actual value),  
Timer: up to 999 hours,  
Heating elements: four-sided heating with large surface area,  
Included:  
Power cord,  
1 pcs. - inner tray with dimensions: 55 x 36 cm,  
Dimensions:  
External: 71 x 55 x 76 cm,  
Internal: 56 x 40 x 48 cm,  
Weight: 50 KG,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

