



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ALC PK110 laboratory centrifuge (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Italian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 150 W,
ALC PK110

is a table top laboratory centrifuge. Its design ensures ease of use, safe operation and a wide range of applications in medical, biochemical, industrial and other laboratories where there is a need to separate the test preparation into components of different densities under the action of centrifugal force,

Technical specifications:

Maximum speed: 4,000 rpm (pivoting equipment), 6,300 rpm (fixed angle equipment),
Maximum RCF at the tip: 2,647 xg (pivoting equipment), 3,661 xg (fixed-angle equipment),
Maximum allowable capacity: 4x135 ml (pivoting equipment), 6x50 ml conical (fixed angle equipment),
Adjustable timer from 1 min to 99 min + continuous time,

Features:

Brushless induction motor,
Microprocessor control,
Digital set/speed readout,
High and low acceleration,
High and low deceleration,
Low noise level,

Included:

User manual in English (PDF),
Power cord,
ALC O-G26 rotor - no.11210021,
4x ALC O-G26/1 cup - no.11210025,
2x ALC adapter no.11172222 (1x 15 ml.),
2x ALC adapter (1x 85 ml.),
Dimensions: 37 x 46 x 28 cm,
Weight: 28 KG,

Has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

