



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EMED apparatus cart 080-040 (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

EMED 080-040

Is a stable and sturdy apparatus cart with free movement,

It has:

Built-in opening door (lockable),

2x cylinder mounts,

Upper top with dimensions: 41 x 40 cm,

Lower shelf with dimensions: 33 x 41 cm,

Dedicated air hoses for connecting cylinders with apparatus,

Includes:

FARUM O2 reducer,

Dedicated key for opening the door,

Overall dimensions: 41 x 48 x 96 cm,

Total weight: 25 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

