



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STATIM 5000 G4 Wi-Fi Autoclave (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Canadian production, (2018),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1300 W,
Autoclave STATIM 5000 G4 ENDO
sensationally performs working in extreme conditions in dozens of hospital operating theaters. The offered model has a built-in WiFi modem, allowing the autoclave to be seamlessly connected to the existing computer network,
The device is equipped with a color touch screen. Such a solution is very user-friendly. It allows quick start of the sterilization cycle, its ongoing control and very easy management of the extensive menu,
Technical specifications:
Class B,
Operating volume - 57 dB,
Cycle recording - built-in Data Logger,
Cycle for packaged 134°C - 18 min,
Cycle for unwrapped 134°C - 9 min,
Plastics and rubber 121°C - 24 min,
Drying - minimum 60 min,
Sterilization and drying of packaged - about 82 min,
Sterilization and drying of unpackaged - about 45 min,
Sterilization of batches of type - A and B,
Cartridge capacity - 5.1 l,
Water tank capacity - 4 l,
LAN operation : WiFi or RJ45 network cable,
Equipped with RS232 port,
Pre-vacuum: thermodynamic,
Drying: thermal,
Features:
Color LCD touch screen display,
Advanced G4 operating system,
Easy operation,
Built-in cycle recorder Data Logger,
Recording on USB-connected flash drive,
LAN network interface,
Energy efficiency,
Fast sterilization cycles - up to 9 min for unwrapped instruments,
Low consumption of distilled water,
Simple design for easy service and periodic maintenance,
Compliance with EN13060 standard for sterilization of type A and B loads,
Selectable menu in the following languages: PL, ENG, CS, DUN, ESP, FR, HUN, IT, KOR, NL, PT, RUM, RUS, UA,
Can be retrofitted with either STATIM cassette or STATIM ENDO cassette at additional cost (cassettes shown in the last photo),
Included:
User manual in Polish (PDF),
Power cord,
16GB USB memory stick,
Bottle with condensate hose,
Overall dimensions: 58 x 41 x 20 cm,
Weight: 40 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

