



## **FAMED ŻYWIEC rehabilitation bed LR-03.9 (Ex-Display)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Polish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

FAMED LR-03

.9 is a comfortable rehabilitation bed with electric height adjustment,

Technical specifications:

Electric adjustment by wired foot control:

Bed height adjustment in the range: 38 - 78 cm,

Back segment inclination angle in range: 70°,

Flexion of the two-segment leg rest in the range: 26°,

Features:

Sturdy and stable construction for safety and comfort,

Folding handrails on both sides of the bed,

Bumpers in each corner,

Dimensions:

Backrest: 78 x 70 cm,

Leg rest segment: 76 x 90 cm,

Overall: 100 x 210 x 73 / 113 cm,

Weight: 120 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

