



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CHEIRON VICTORIA II VERSA surgical suction machine (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- Czech production,
- Power supply: 230V,
- Frequency: 50/60 Hz,
- Power: 110VA,
- CHEIRON VICTORIA II VERSA is a high-performance model. Versatile, designed especially for wards, operating rooms, postoperative rooms, intensive care units, obstetrics and endoscopic applications,
- Specification:
 - Maintenance-free unit designed for continuous operation,
 - The casing is made of very durable, shatterproof, fireproof material with UV protection,
 - The pump body is made of aluminum alloy,
 - Overfill protection,
 - Safe vacuum regulator,
 - Smooth suction power regulation,
 - Illuminated on/off switch,
 - Capacity: 38 l/min,
 - Max negative pressure: 90 kpa,
 - Drainage: 0 - 40 cm H₂O,
 - Noise and vibration dampening system for exceptionally quiet and smooth operation at an incredibly low noise level of 39.4 dB
- Included in the set:
 - Reusable 1 L jar from SERRES (autoclavable 121°C) – has a new disposable SERRES 1 L insert,
 - Power cord,
- Dimensions: 42 x 26 x 34 cm,
- Weight: 12 KG,
- Cat 02,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
 - 6 months for domestic market (Poland),
 - 3 months to international market,
 - Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
- Financing options (Poland only): Installments, Leasing, Loan,

