



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ADEC 1040 dental unit (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Upholstery will be replaced with a new one in the color chosen by the customer from our color chart,

Dental unit
ADEC 1040

is a compact, ergonomic dental unit with electromechanical drive,
Armchair operator's table placed on pantograph arm:

5 upper silicone sleeves:
Water-air blower (three-way),
Micromotor (with light)
NSK Ti205L
(speed: 18,000 - 22,000 RPM),
NSK Ti-MAX X600L Turbine:
Optical fiber made of cellular glass,
Ceramic bearings,
Clean head system,
Microfilter,
Push-button drill bit replacement,
For quick release
LED (NSK PTL-CL-LED)

,
Handle for scaler
WOODPECKER
(autoclavable 135°) - has 4 tips and a wrench,
Polymerization lamp
TOP LIGHT

,
Smooth control of the operation of pneumatic tips with a foot control,
Free manual extension of the chair operator's auxiliary tray,
Operator lamp
MOCOM 12 V / 50 W
with free vertical and horizontal adjustment,
Electric seat adjustment by foot control:
Seat height adjustment in the range: 48 - 75 cm,
Backrest tilt: -15° / 0° / 70° (automatic seat back synchronization),
4 pre-programmed patient functions on foot control,
Possibility to unfold the chair to Trendelenburg position,
Freely manually adjustable headrest height and tilt,
Includes:

NSK Ti-MAX X600L Turbine,
NSK Ti205L micromotor (with light),
NSK PTL-CL-LED Quick Release,
WOODPECKER scaler handle(4 tips and wrench),
TOP LIGHT polymerization lamp,
Water and air blower (three-way),

Dimensions:
Headrest: 20 x 18 cm,
Backrest: 65 x 45 cm,
Seat: 58 x 120 cm,
Total: 90 x 130 x 194 cm,
Weight: 170 KG,

Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:

6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,
The warranty period does not include bulbs and other light sources installed in the device,
Possibility of financing (Poland only): Installments, Leasing, Loan,

