



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ECG device WELCH ALLYN CP 200 electrocardiograph (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230V,
- Emergency power supply: via 4 new built-in batteries 6.4 V / 6 Ah,
- Frequency: 50/60 Hz,
- **The Welch Allyn CP 200** is a 12-lead resting ECG system that can also perform spirometry testing. The CP 200 is equipped with an alphanumeric keyboard and an LCD display. The display can show a quick color preview of the ECG reading along with the ability to edit settings. The Welch Allyn system provides full-page (8.25 x 11) ECG printouts and can store up to 50 ECG and 50 spirometry measurements in its internal memory. The ECG system offers the ability to transfer patient data via SD card, UBS system or your facility's EMR system,
 - High-resolution thermal printer generates full-size standard and custom reports,
 - Full alphanumeric keyboard and LCD display for quick and easy entry of patient information,
 - Advanced filters for optimal ECG recording quality ensure accurate readings,
 - Configurable report formats allow you to print information the way you want,
 - Battery power allows easy transport of the device to the patient,
 - The instant power-on function allows you to quickly turn on the device and immediately perform an ECG test,
- **Technical specifications:**
 - Resting 12-lead ECG,
 - LCD screen for easy viewing,
 - Performs spirometry tests,
 - Stores up to 50 ECG readings and 50 spirometry readings,
 - Full-page prints (8.25 x 11),
 - Configurable report formats,
- **Possibility to select menu language: EN, DE, FR, ES, IT, PT, SE, NR, NL, FI, PL, RU, GR, HR, ML,**
- **Control panel keyboard in English,**
- **Included in the set:**
 - English user manual (PDF),
 - Power cord,
 - Mobile transport trolley,
 - 10-lead ECG cable,
- Dimensions:
 - Overall: 41 x 40 x 16 cm,
- Weight: 10 kg,
- Cat 03,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

