



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical monitor KARL STORZ 20090331 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Screen: 15" inches,
Touch LCD display: 1024 x 768 px,
Inputs:
DVI(IN),
VGA,
RS-232,
Selectable menus in EN, FR, DE, IT, ES, RU, PT,
Included:
Power cord,
Power supply,
Dimensions:
Overall: 37 x 30 x 6 cm,
Weight: 2 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

