



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



GE DASH 5000 Patient Monitor/Cardiac Monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V,
Emergency power supply: through two built-in batteries 11.1 V / 3.5 Ah
(not included, can be retrofitted with new ones at additional cost),
Frequency: 50/60 Hz,
Power: 75 W,
Specifications:
12.1" color LCD display with resolution: 800 x 600 px., Provides an efficient screen image from a distance,
Slots: Defib Sync. / LAN / Aux,
Custom display options
enable:
View, and analyze data in graphical and numerical format,
Compilation of real-time measurements and trends,
Free configuration of each screen element - from curves to parameter names,
Fully flexible alarm settings for personalized monitoring,
Bright bedside alarm light helps locate the patient immediately in critical situations,
Both baseline and side stream CO measurement options.

2

low-flow for applications ranging from neonatal to adult patients,
Lightweight, durable and ergonomic design, optimized for demanding transport,
Features:
TRIM KNOB navigation dial operation for easy information entry and screen navigation,
Easy-to-use interface,
Equipped with:
NIBP - New NIBP cuff with dual air tubing,
ECG - New 3-lead ECG lead,
SOP

2

- New clip sensor,
TEMP - New temperature sensor
(Not included - can be retrofitted at additional cost)

,

PI1/3 and PI2/4
- Measurement line can be retrofitted at additional cost,
Selectable menu in EN, FR, IT, ES, DE, NL, SE, DK, HU, NR, PA, PL, RU, CZ, JP,
Included:
User manual in English (PDF),
Power cord,
NIBP - New NIBP cuff with dual air leads,
ECG - New 3-lead ECG cable,
SOP

2

- New clip sensor,
Dimensions: 22.5 x 31.5 x 30 cm,
Weight: 10 KG,
Cat 01,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

