



PRAISTON Sp. z o.o.
ul. Górska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



POPPELREUTER Test Board (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Poppeleruter board
is used to assess concentration, metastability and divisibility of attention under conditions of imposed time pressure,
Assessment of concentration, divisibility and metastability of attention,
Rigid boards for hanging or standing,
Durable and clear printing,
Includes:
Mounting stand,
Board dimensions: 65 x 65 cm,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

