



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NIDEK NT-2000 non-contact tonometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 90 VA,
Technical specifications:
Measurement range: 1-60 mmHg,
Built-in 5" CRT display,
Working distance of 11 mm (from the air nozzle to the cornea of the eye),
Automatic pupil guidance,
Built-in thermal printer (paper width 56 mm),
Electric canopy lift,
Automatic safety interlock to prevent contact with the patient's eye,
Computer interface for data transmission,
Manual chin adjustment within 7 cm,
Menu in language: English,
Included:
User manual in English (PDF),
Power cord,
Dimensions:
Overall: 26 x 52 x 45 cm,
Weight: 19 kg,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,



