



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FAMED ŻYWIEC WS-08.5 support roller (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Polish production,
- **The WS-08.5 support roller** is used to support the patient in a lateral position,
- Construction made of stainless steel,
- Possibility of adjusting the rotation of the roller in the range of 360°,
- **Included in the set:**
 - Bone attaching to the rail,
- Dimensions:
 - Roller: 29 x 30 cm,
- Weight: 2.5 KG
- Cat 03,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

