



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FUKUDA DENISHI CARDIMAX FX-7402 ECG machine (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Emergency power supply: via built-in new 3800 mAh / 9.6 V battery,
Frequency: 50/60 Hz,
Power: 80 VA,
The apparatus is ideal for private practices and doctors' offices, equipped with LCD screen for easy operation and membrane keys,
Features:
Clear, touch-screen ECG display,
Built-in memory and 12 leads with real-time ECG waveforms,
A4 recording format,
Built-in thermal printer,
Technical specifications:
LCD display 320 x 240 px,
ECG signal: 12 standard leads,
Test sensitivity: $\frac{1}{4}$, $\frac{1}{2}$, 1.2 cm/mV - or automatic,
Recording speed: 5 / 10 / 25 / 50 mm/sec,
Digital filters: 0.05 to 150 Hz,
Measurement operations:
Heart Rate / R-R time / P-R time / QRS time / QT time / Qtc / Electrical Axis / SV1 / RV5(6),
Heart rate measurement: from 20 to 300 bpm (accuracy +/- 2 bpm),
Selectable menu in ENG, GER, FR, IT, ESP,
Instrumentation placed on a mobile cart:
Shelf designed for accessories,
Built-in 3 outlets for powering accessories,
Included:
User manual in English (PDF),
Power cord,
Dedicated cart,
10-lead ECG cable,
Overall dimensions: 45 x 40 x 95 cm,
Weight: 30 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

