



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDINSTRUS Ophthalmic Chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Lithuanian production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 50 W,
Upholstery with medical certificate in brown color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),
Armchair
MEDINDUSTRIES
is a multifunctional chair ideally suited for ophthalmology offices,
Technical specifications:
Electric height adjustment of the chair seat in the range of 55 - 70 cm,
Manual adjustment of the chair backrest with a lever in the range of 75°,
The chair has the ability to rotate around its own axis in the range of 360° (with the possibility of locking in the selected position),
Retractable patient armrests on both sides,
Retractable patient step,
Included:
Power cord,
Wired foot control,
Dimensions:
Backrest: 55 x 75,
Seat: 55 x 48 cm,
Overall: 70 x 74 x 90 cm,
Weight: 75 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

