



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDINTRUS cosmetic and treatment chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Lithuanian production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 50 W,  
Upholstery with medical certificate in brown color (it is possible to replace the upholstery with a new one, from our color pattern for an additional fee),  
Armchair  
MEDINDUSTRIES  
is a multifunctional chair ideally suited for beauty and treatment rooms,  
Technical specifications:  
Electric height adjustment of the seat of the chair in the range of 55 - 70 cm,  
Manual adjustment of the backrest of the chair with a lever in the range of 75°,  
The chair has the ability to rotate around its axis in the range of 360° (with the possibility of locking in a selected position),  
Retractable patient armrests on both sides,  
Retractable patient step,  
Included:  
Power cord,  
Wired foot control,  
Dimensions:  
Backrest: 55 x 75,  
Seat: 55 x 48 cm,  
Overall: 70 x 74 x 90 cm,  
Weight: 75 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

