



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NVACARE PLATINUM S oxygen concentrator (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

P.R.C. production,

Power supply: 230 V,

Features:

Increasing blood oxygen levels,

Reducing the progression of respiratory diseases,

Improving well-being,

Increasing physical fitness,

Improving oxygenation of internal organs,

Increasing the body's immunity,

Improving concentration,

Alleviation of ailments: weakness, sleep disorders, inflammation, headaches and dizziness,

Specifications:

Oxygen flow control (guarantees a very high level of oxygen concentration in the full flow range from 0.5 l to 5 l.),

Oxygen concentration of 87% ,

Easy to transport (light weight, handle and 4 maneuverable wheels improve daily comfort),

Easy to use (it has easy access to the humidifier and air filter, and intuitive adjustment of oxygen flow with a dial. To make it work, all you need to do is plug it in and start it up),

Includes:

Power cord,

Overall dimensions: 38 x 30 x 68 cm,

Weight: 21 KG,

Has a current inspection and is ready to work,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,



