



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



LEC Laboratory Refrigerator LM-PGRC273UK (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
EU production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Professional laboratory refrigerator
LEC LM-PGRC273UK
with a capacity of 273 liters is designed from the ground up for safe storage of drugs in the cold chain,
Specifications:
Temperature control range: +2 °C to +8 °C,
Built-in temperature alarm,
Internal refrigerator lighting,
Automatic defrost system,
Controller temperature measurement accuracy up to 0.1°C,
Antimicrobial handle helps prevent bacterial growth and minimizes the risk of cross-contamination,
High-Tech controller shows the temperature of the air and simulated load, allowing you to check the temperature of the contents stored in most models,
Control panel lock to prevent accidental changes to selected settings,
Open door sensors warn when doors are left open,
Advanced controller to display min and max temperature and manage warning alarm system,
Internal light , which can be turned on and off using an external control panel,
Refrigerant: R600a,
Included:
5 pcs. - shelf (36 x 40 cm),
Dimensions:
External: 60 x 62.5 x 156.6 cm,
Internal: 48 x 44.5 x 138.2 cm,
Weight: 107 KG,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

