



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **LEC LM-PGRC273UK Laboratory Refrigerator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- EU production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The professional LEC LM-PGRC273UK** laboratory refrigerator with a capacity of 273 liters is designed from the ground up for the safe storage of medicines in the cold chain,
- **Specification:**
  - Operating temperature range: +2 °C / +8 °C,
  - Built-in temperature alarm,
  - Internal refrigerator lighting,
  - Automatic defrosting system,
  - Controller temperature measurement accuracy up to 0.1°C,
  - Antibacterial handle helps prevent bacterial growth and minimizes the risk of cross-contamination,
  - The High-Tech controller displays the air and simulated load temperature, allowing you to check the temperature of the stored contents in most models,
  - Control panel lock to prevent accidental changes to selected settings,
  - Door open sensors warn when a door is left open,
  - The advanced controller allows you to display min and max temperatures and manage the warning alarm system,
  - Internal light that can be turned on and off via an external control panel,
  - Coolant: R600a,
- **Included in the set:**
  - 5 pcs. – shelf (36 x 40 cm),
- Dimensions:
  - External: 60 x 62.5 x 156.6 cm,
  - Internal: 48 x 44.5 x 138.2 cm,
- Weight: 107 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

