



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEC PG-907 laboratory refrigerator (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

EU production,

Power supply: 230 V,

Frequency: 50/60 Hz,

LEC PG-907 professional laboratory refrigerator with a capacity of 273 liters is designed from the ground up for safe storage of drugs in the cold chain,

Specifications:

Temperature control range: +2 °C to +8 °C,

Built-in temperature alarm,

Internal refrigerator lighting,

Automatic defrost system,

Refrigerant: R600a,

AUX connector,

Included:

1 pc. - door key,

5 pcs. - shelf (36 x 40 cm),

Dimensions:

External: 60 x 62.5 x 156.6 cm,

Internal: 48 x 44.5 x 138.2 cm,

Weight: 107 KG,

It has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

