



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCOPEASSIST PFE-901 irrigation pump (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
U.K. production, (2019),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 360 VA,
SCOPEASSIST PFE-901
is a powerful flushing pump that guarantees uninterrupted irrigation with pressure control,
Constant patient monitoring,
Thoughtful design of the roller element ensures easy installation of the tube with integrated control sensor,
Fast start-up time,
Enhanced safety for pressure and flow,
Digital control panel,
3-stage flow control,
The pump has a dedicated wall mount ,
Includes:
Foot attachment
Power cord,
Dedicated wall mount,
Dimensions: 40 x 25 x 12 cm,
Weight: 6 KG,
Has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

