



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Video-laparoscope OLYMPUS HD II ENDOEYE WA50042A (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German manufacture,
Olympus Endoeye HD II
is a video - laparoscope with 1080i High Definition Television (HDTV) display standard,
With 1080 scan lines, the image resolution is about twice that of a conventional TV display,
Capable of working with Narrow Band Imaging (NBI) technology, which can be used to support diagnostics,
Specially designed, distally mounted 1080i H-CCD-compatible image sensor provides images with clarity and resolution that enable even the smallest details to be seen,
Advanced optical design provides greater depth of field, eliminating the need for manual focusing and providing consistently excellent images with a sharp center and sharp edges in 16:9 widescreen format,
The anti-fog function provides a clear image throughout the procedure,
Durable and ergonomic autoclavable design,
One-handed controllability with three programmable remote control buttons,
Specifications:
Viewing angle: 30°,
Outer diameter of endoscopic probe tip: 10 mm,
Working length: 33 cm,
Total length: 500 cm,
Included:
OLYMPUS WA50042A autoclavable sterilization container,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

