



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KNÜRR Endoscope Trolley (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Specifications:
Top drawer lockable,
Power switch on the front of the cart,
Equipped with electrical installation,
Endoscope holder with height adjustment in range: 60 cm,
Removable rear housing,
Mobile on 4 casters including 2 with brake,
Included:
Power cords (4 pcs),
Dimensions:
Top: 42 x 44 cm,
Shelves: 45 x 45 cm,
Top drawer: 34 x 37cm,
Lower drawer: 36 x 38cm,
Overall dimensions: 60 x 55 x 122 cm,
Weight: 49kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

