



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



KNÜRR Endoscope Trolley (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

German production,

Specifications:

Top drawer lockable,

Power switch on the front of the cart,

Equipped with electrical installation,

Endoscope holder with height adjustment in range: 60 cm,

Removable rear housing,

Mobile on 4 casters including 2 with brake,

Included:

Power cords (4 pcs),

Dimensions:

Top: 42 x 44 cm,

Shelves: 45 x 45 cm,

Top drawer: 34 x 37cm,

Lower drawer: 36 x 38cm,

Overall dimensions: 60 x 55 x 122 cm,

Weight: 49kg,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

