



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DATEX-OHMEDA TEC 7 inhalant vaporizer (Isoflurane/Isoflurane) (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

Manufactured by U.S.A.,

DATEX-OHMEDA TEC 7 vaporizers

designed for the delivery of the anesthetic agent isoflurane, combines significant performance, convenience and ergonomic benefits with important safety features that help improve agent delivery and reliability while reducing overall operating costs,

Provides consistent initial vapor concentration over the entire range of flow rates used in clinical practice,

Easy-to-turn dial and fine adjustment capability helps to precisely control the concentration of anesthetic administered,

The vaporizer dial for isoflurane allows step changes in 0.2% increments for concentrations up to 1% and 1% increments for concentrations from 1% to 8% ,

To reduce the frequency of filling, the TEC 7 vaporizer holds 225 ml of anesthetic.

The TEC 7 vaporizer can be installed in Datex-Ohmeda anesthesia machines equipped with the Selectatec attachment manifold. No tools are required for this operation, and the anesthesia machine does not need to be taken out of service,

TEC 7 vaporizers can be installed on a single mounting manifold with Datex-Ohmeda TEC 4, TEC 5, TEC 6 and TEC 6 Plus vaporizers. The listed vaporizers are compatible with the co-locking mechanism that the TEC 7 vaporizer is also equipped with,

Ability to read the level of liquid anesthetic thanks to an indicator on the front of the evaporator,

Selectatec (S - 2000) attachment is used in the given vaporizer,

Filler type: Easy-Fill,

The vaporizer has signs of use that do not affect the comfort of use,

Included:

User manual in English (PDF),

Power cord,

Overall dimensions: 22 x 11 x 25 cm,

Weight: 8 kg,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

