



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Nasopharyngoscope ALTHEA (3.4 mm) (New)

New,
Realistic pictures of the product,
German manufacture,
Nasopharyngoscope ALTHEA
is characterized by an excellent, clear image and excellent ergonomics,
The tip of the endoscope is tiltable and moves in two planes (within 180°), making the endoscope very handy during the procedure,
Thanks to the use of modern technology, an unprecedented resolution of 18,000 pixels has been achieved,
Technical specifications:
Working length: 300 mm,
Diameter: 3.4 mm,
Field of view: 85°,
Depth of field: 1 - 50 mm,
Tip bend: up: 150° / down: 150°,
Included:
Instruction manual,
Pump for leak test,
Original carrying case for transporting the instrumentation,
Has up-to-date inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 24 months,
Possibility of financing (Poland only): Installments, Leasing, Loan,

