



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



YESON YS-22L-E Autoclave (Ex-Display)

New (post-show),
Real photos of the product,
Manufactured by P.R.C., (2023),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 1950 W,
Modern autoclave
Yeson E series (LCD)
With a capacity of 22 liters,
European class B in accordance with the EN 13060 standard,
The device is designed for use in: hospitals, medical facilities (in clinics), including medical and dental offices, as well as laboratories. The autoclave can also be used for non-medical and veterinary applications,
Features and advantages of the E-series autoclave:
Stainless steel chamber,
Class B fractional vacuum,
Built-in printer,
USB connector for flash drive,
LCD display,
Touch control buttons,
Polish and English menus,
Lock against accidental opening,
Automatic fuses,
9 operating programs,
2 test programs,
Delayed start function,
Pre-heat function,
Password protection function,
Technical specifications:
Medical device: class IIb,
Chamber capacity: 22 liters,
Internal tanks for distilled and used water,
Loading method: basket with 3 trays,
Power supply: 230V \pm 10% , 50/60 Hz,
Rated power: 1950 W,
Maximum working pressure: 0.21 ~ 0.24 MPa / 2.1 bar ~ 2.4 bar,
Maximum operating temperature: 134°C ~ 137°C,
Included:
3x tray of 18.5 x 39 cm,
Tray holder,
Drain hose,
Power cord,
Dimensions:
Dimensions with door closed: 69.5 x 50.5 x 42.5 cm,
Chamber dimensions: 25 x 45 cm,
3x tray dimensions: 18.5 x 39 cm,
Weight: 60 KG,
It has a current inspection and is ready to work,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

