



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WINNCARE HELIOS 150 EU patient lift (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Danish production, (2021),  
Power supply: 230 V,  
Emergency power supply: via 2 new built-in 12 V / 2.9 Ah batteries,  
Frequency: 50/60 Hz,  
Helios 150  
is a modern mobile lift with a compact design, ideal for use in home care, nursing homes and hospitals,  
With its lightweight design and ability to be easily disassembled into two parts, the Helios 150 provides exceptional mobility and convenience of use, without the need for tools,  
Helios 150  
is designed for convenient and safe patient lifting, and its mobility also allows for short transport. It is an indispensable tool in the daily work of caregivers, nurses and medical personnel,  
Features:  
Ergonomics and safety: The design ensures that the patient's center of gravity is above the low-friction wheels, which greatly facilitates maneuverability and ensures safety during lifting,  
Functionality: The lift is equipped with a manually extendable running base, which allows for quick folding and unfolding, adapting the device to different conditions of use,  
Efficiency: The capacious battery ensures long-lasting operation and allows multiple lifts to be carried out on a single charge,  
Reliability: In emergency situations, the lift has both electric and manual lowering options for safer operation,  
Specifications:  
Lifting range min/max: 42 / 164 cm,  
Lifting method: electric,  
Leg unfolding method: manual,  
Turning radius: 130 cm,  
Front/rear wheel diameter: 75/100,  
Included:  
Power cord,  
Universal transport suspension size M,  
Dimensions:  
Lifting compartment min/max: 42 / 164 cm,  
Base outer width min / max: 61 / 95.6 cm,  
Base inner width min / max: 39 / 51 cm,  
Mounting base height min / max: 10 / 15 cm,  
Maximum load capacity: 150 KG,  
Weight: 35 KG,  
Has an up-to-date inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

