



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Arm support (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real product photos,
- Possibility of free adjustment of support,
- Fastening belt,
- Dimensions: 11.5 x 40 cm,
- Weight: 2 KG,
- It has a current inspection and is ready to work,
- Issued Technical Passport (Service Report) valid for 12 months,
- Guarantee:
 - 6 months for domestic market (Poland),
 - 3 months to international market,
 - Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,
- Financing options (Poland only): Installments, Leasing, Loan,

