



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NIKON NS-1 Slit Lamp (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Light source: halogen bulb: 12V / 30W,
Technical specifications:
Infinitely adjustable light intensity,
Magnification eyepiece: 15 x/14,
Infinitely adjustable magnification change from 7.5 to 30,
Infinitely adjustable focus,
Stepless adjustment of slit width: 0-20 (by scale),
Stepwise adjustment of slit length: 0.2 / 1 / 2 / 5 / 10 / 14 + blue filter (according to the scale),
Stepless adjustment of slot length: 1 - 12,
Slit angle adjustment: 0 - 180°,
Infinite adjustment of the patient's chin in the range: 0 - 7 cm,
Vertical tilt angle of the lamp (5° gradual adjustment),
Rotatable eyepiece (lamps): 90° right / left,
Infinitely adjustable up and down lamp (joystick) in the range: 0-3 cm,
Fixation point,
Built-in 5 types of filters:
Transparent filter,
Heat absorption filter,
ND filter (28% transmission),
Green filter,
Blue filter,
The device placed on a mobile KEELER table with manual height adjustment from 65 - 88 cm,
Included:
Power cord,
Power supply,
Dimensions:
Slit lamp: 33 x 40.5 x 70 cm,
Overall: 60 x 45 x 130/153 cm,
Weight: 14.5 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

