



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



ASTAR ETIUS electrotherapy device (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production, (2013),
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 70 W,
ETIUS electrotherapy apparatus
has two fully independent electrotherapy channels with the possibility of coupling for four-electrode treatments. The output circuits can operate in two modes: patient current stabilization or output voltage stabilization. The latter is used for non-stationary and combined techniques. The device has a built-in database of ready-made therapeutic procedures, and there is also the possibility for the user to create his own treatment programs

The device offers completely new possibilities: the application of microcurrents, unipolar wave currents and specially developed TENS currents for the treatment of spastic paralysis,

Specially designed graphical menu makes operation of the device easy and intuitive,

Features:

Two completely independent treatment channels,

Operation mode: program/manual,

Large, easy-to-read display,

Graphical mode operation,

Database of built-in treatment programs,

Selection of disease units by name,

Treatment timer,

Electrode test,

Statistics of performed treatments,

Adjustment of the volume of the sounder,

Technical specifications:

Interferential, Kotz unipolar waveform: 100 mA,

TENS: 140 mA,

Diadynamic / pulsed: 60 mA,

Galvanic: 40 mA,

Tonolysis: 100 mA,

Microcurrent: 1000 μ A,

Maximum patient circuit voltage amplitude (CV mode): 140 V,

No electrodes included - electrodes can be retrofitted at additional cost,

Instrumentation placed on a dedicated ASTAR mobile cart,

Includes:

User manual in Polish (PDF),

Overall dimensions: 65 x 55 x 100 cm,

Total weight: 25 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

