



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WELCH ALLYN PIC50 Defibrillator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Manufactured by U.S.A.,
Power supply: 230 V / DC 12 V (car socket),
Emergency power supply: through a new built-in battery,
Frequency: 50/60 Hz,
Device
Welch Allyn PIC 50
is a reliable vital function monitor, defibrillator and pacemaker,
Intelligently designed interface and controls make it easy to perform rescue operations,
Technical specifications:
Defibrillator:
Biphasic,
Energy selection:
External: 2, 5, 7, 10, 20, 30, 50, 70, 100, 150, 200, 300, 360J,
Internal: 2, 3, 4, 5, 6, 7, 8, 9, 10, 15, 20, 30, 50J,
Charging time:
Charge indicator: audible and graphic,
(Portable) has a sturdy handle to carry the defibrillator freely,
Features:
SYNC - Synchronization (cardioversion),
PACER - Pacemaker,
Built-in display:
Size: diagonal 6.5 inches (16.5 cm), no glare,
Type: active matrix color TFT display,
Display resolution: 640 x 480 pixels,
Built-in thermal printer,
Paper width: 50 mm,
Print speed: 25 mm/s, 50 mm/s,
Menu in language: Polish,
Keyboard on the control panel in: English,
Equipped with the following modules and cabling:
ECG - 4-lead ECG cable,
SPO2 - No sensor
(can be retrofitted for an additional charge),
NIBP - Cuff (23-33 cm) + air line,
CO2 - No measurement line
(can be retrofitted for an additional charge),
TEMP - No sensor
(can be retrofitted for an additional charge),
IBP x 2 - Invasive pressure measurement - no transducer cable
(can be retrofitted for an additional charge),
Included:
User manual in English (PDF),
Power cord,
Disposable multifunctional EURODEFIPADS electrodes,
Dedicated carrying case,
Dimensions:
Total instrumentation: 26 x 35 x 37 cm,
Weight:
Defibrillator: 4.5 KG,
Total Instrumentation: 6 KG,
Cat D145,
It has a current inspection and is ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

