



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STOLTER FG-01 gynecological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Polish production,
Gynecological chair FG-01
- is designed for gynecological and urological examinations and procedures, in health care facilities,
Construction made of steel profiles covered with powder varnish,
New upholstery with medical certificate in blue color (it is possible to replace the upholstery with another one from our color pattern for an additional fee),
Specification:
Backrest adjustment in the range: 0 - 60°,
Chair tilt:
Trendelenburg approx. 0 - 25°,
Anti-Trendelenburg approx. 0 - 10°,
Possibility to obtain the reclining position,
Fixed seat height of the chair,
Tilt of the chair in the reclining position: 0 / - 30°,
Equipped with:
Stainless steel rails for attaching accessories,
2x removable patient handle,
2x removable bolster,
2x leg rest,
Pull-out stainless steel bowl with dimensions: 32 x 26 x 10cm,
Paper holder with a width of 33cm,
In the base of the chair folding platform,
Features:
Adjustable knee pads with large range,
Possibility to change the position of the patient holder,
Dimensions:
Backrest: 62 x 85 x 6 cm,
Seat: 62 x 37 x 6 cm,
Base: 60 x 100 cm,
Overall: 146 x 80 x 158 cm,
Weight: 75.5 kg,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

