



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MITSUBISHI P66E Videoprinter (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
Japanese production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Videoprinter Mitsubishi P66E
is designed to print images from video devices including medical diagnostic equipment, CCTV cameras and VCRs with full photographic quality,
High-speed, high-resolution thermal printer for medical applications,
Specifications
:
Built-in display,
Retractable control panel,
Monochrome printing,
Scaling : 256 levels,
Resolution : horizontally 1216 pixels x vertically 600 lines PAL (standard),
Print speed: 5.8 sec,
Video in/out,
Included:
User manual in English (PDF),
Power cord,
Dimensions: 15.4 x 31.4 x 12.5 cm,
Weight: 3.9 KG,
It has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for a period of 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

