



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS CV 145 video processor with CLE-145 light source (set) (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Japanese production,  
The set includes:  
OLYMPUS Model CV-145 EVIS EXERA image processor:  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Automatic white balance control,  
It has the following inputs and outputs:  
MONITOR REMOTE,  
PRINTER REMOTE, D. FILE REMOTE,  
USB 2.0 - 2 pcs for keyboard connection,  
MONITOR OUT,  
COMP OUT,  
RGB OUT,  
Y/C OUT,  
LIGHT CONTROL,  
Features FREEZE function,  
Electronic edge sharpening of endoscopic images improves imaging and diagnostics,  
Adjustable color saturation: red and blue,  
Foiled control panel,  
Dimensions: 46 x 36.5 x 7.5 cm,  
Weight: 8 kg,  
OLYMPUS Model CLE-145 EVIS EXERA light source:  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Adapted for flexible endoscopes,  
Main and emergency lighting: 2 150 W halogen bulbs,  
Easy access to the bulb from the side of the unit,  
Lamp life: approx. 50 h with continuous use,  
Aperture control capability,  
Features LIGHT CONTROL port,  
Foiled control panel,  
Automatic or manual adjustment of image brightness,  
Built-in air pump with electronic power control,  
Dimensions: 45 x 36.5 x 15 cm,  
Weight: 13 kg,  
OLYMPUS MAJ - 844 keyboard,  
Light Control cable MH - 966,  
Bottle for endoscope flushing,  
Pigtail cable MAJ - 843,  
Power cords - 2 pcs,  
Overall dimensions: 45 x 36.5 x 22 cm,  
Weight of the set: 21 kg,  
Version B  
It has a current Technical Passport issued,  
Warranty:  
6 months for domestic market - Poland,  
3 months for the international market,  
Possibility of financing (Poland only): Installments, Leasing, Loan,

