



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DERUNGS HALUX N50 treatment lamp (Ex-Display)**

New (post-show),  
Real photos of the product,  
German manufacture, (2020),  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 13 W,  
Lamp  
Halux N50  
is ideal for small procedures and general examination,  
Technical specifications:  
Fast cleaning due to closed design,  
Maintenance-free - High-quality materials and durable LED technology,  
Safe and easy positioning,  
High performance - 50,000 lux from 0.5 m,  
Large, uniform light field for maximum viewing comfort,  
Treatment field - Ø 18 cm from 0.5 m,  
Color temperature Fixed: 4 400 K,  
Articulated arm with free adjustment,  
Integrated power supply,  
Included:  
Side rail mount,  
Dimensions:  
Total length of articulating arm: 85 cm,  
Weight: 1.3 KG,  
Has a current inspection and is ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty at a surcharge to 12-months for domestic and international markets,  
The warranty period does not include bulbs and other light sources installed in the device,  
Possibility of financing (Poland only): Installments, Leasing, Loan,



