



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH footrests (knee rests) type Goepel (set) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- They are used to support the lower limbs in a lying or sitting position,
- They are equipped with polyurethane mattresses with dimensions: 16 x 27.5 cm,
- The footrests mounted on a ball joint allow for adjustment in many planes,
- The set includes two footrests with arms and two mounting bones for the rail with adjustable height and angle,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for domestic market (Poland),**
 - **3 months to international market,**
 - **Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,**
- Financing options (Poland only): Installments, Leasing, Loan,

