



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **KEELER PULSAIR DESKTOP non-contact tonometer (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real photos of the product,

U.K. production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 30 VA,

Pulsair Desktop

is an electric, non-contact screening tonometer for measuring intraocular pressure, which is one of the key tests in diagnosing glaucoma,

The device has a color display that facilitates simplicity and improves user comfort. After calibrating the system, the tonometer will automatically take IOP measurements. The system also allows the user to print IOP measurements,

Specifications:

Calibration range: 5 mmHg to 50 mmHg,

Working range: 20 mm from the lens to the patient's cornea,

Display scale: Single 16-character rows (alphanumeric display),

Illumination system: Infrared LED,

Menu in language: English,

Included:

User manual in English (PDF),

Power cord,

Dimensions: 41 x 24.5 x 47.5 cm,

Weight: 16 KG,

Has a current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for a period of 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

