



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MAQUET RADIUS gynecological chair 1557.01BE (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
German production,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 290 VA,
MAQUET RADIUS
is used for gynecological and urological examinations and procedures,
Upholstery with medical approval in green color (it is possible to replace the upholstery with a new one from our color pattern for an additional fee),
Electromechanical adjustment by foot attachment and wired handheld remote control:
Seat height in the range of 60 - 90 cm,
Tilt angle:
Back segment within (- 9°) / 0° / 57°,
Seat in the range of 18°,
When adjusting the angle of the back segment, the seat angle is automatically adjusted,
3 preset buttons with function memory,
STOP function - button with activation of immediate interruption of the function being performed,
The chair has:
Storage compartment for a 38 cm wide roll of backing paper,
A cutout in the upholstery of the back segment for the backing paper,
Retractable waste bowl measuring 18 x 25 cm,
Reversible footrests,
Ability to unfold the chair to the reclining position,
Includes:
Foot attachment,
Wired remote control,
Waste bowl,
Dimensions:
Back segment: 91 x 61 cm,
Seat: 61 x 31 cm,
Length in extended position: 175 cm,
Overall: 158 x 75 x 138/168 cm,
Maximum load capacity: 180 KG,
Weight: 108 KG,
Has a current inspection and is ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6-months for domestic market (Poland),
3-months for the international market,
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,
Financing options (Poland only): Installments, Leasing, Loan,

