



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CRONOS TP-B1 Portable Oxygen Concentrator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Real photos of the product,
P.R.C. production,
Power supply: 230 V,
Emergency power supply: through 1 built-in new Li-Ion 14.8 V battery,
Frequency: 50/60 Hz,
Power: 32 W,

Oxygen concentrator is the most reliable, effective and convenient oxygen source today. The concentrator separates oxygen from the air and delivers pure oxygen of high concentration to the patient through a tube,

The concentrator is especially useful for people who use oxygen therapy while traveling, and who lead active lifestyles. When driving, just plug in the car charger and the concentrator can run for many hours without interruption,

Features:

Flow rate averages 3 L/min. (+/-1.5 L/min depending on the need),

High oxygen concentration,

Alarm system for power cut-off and too high or low pressure,

Compressor operation protection,

Small size and weight,

Low power consumption,

Included:

Instruction manual,

Power supply,

Battery power supply,

14.8V battery,

Vaporizer tube,

Dimensions: 26 x 11 x 17 cm,

Weight: 2.5 KG,

Has a current inspection and is ready for operation,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6-months for domestic market (Poland),

3-months for the international market,

Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,

Financing options (Poland only): Installments, Leasing, Loan,

